

35 years of experience
in traceable sample
storage




Class 7 clean room
injection molding
and assembly

Production in
the Netherlands and
the United States

MICRONIC TUBE STORAGE RACK RANGE

Roborack-96

Holds 96 Individual Tubes up to 2.50ml!



Superior features

- Based on the Society for Laboratory Automations and Screening (ANSI/SLAS) standards for storage racks
- Standard available with laser etched 1D rack barcode
- Secure and safe locking system
- Twist-lock
- Tapered walls allow precise tube placement into racks (sure shot design)
- Open rack bottom
- Compatible with the Micronic Code Readers
- Class 7 cleanroom production
- RNase/DNase and Pyrogen free

Ordering information

Catalog number	Description	Unit
MPW51001BC3	Roborack-96 white, laser etched (including standard cover) (barcoded A1-H1 side)	10 racks
MPW51012BC3	Roborack-96 white, laser etched (including high cover) (barcoded A1-H1 side)	10 racks
MPW51016BC3	Roborack-96 white, laser etched (including extra high cover) (barcoded A1-H1 side)	6 racks
MPW51017BC3	Roborack-96, white, laser etched (including low cover) (barcoded A1-H1 side)	10 racks



To contact Micronic, email to sales@micronic.com or call us at:

Micronic America (United States, Latin America and Canada): +1 484 480 3372 | **Micronic Asia** (Asia and Pacific): +65 85 028 878 | **Micronic Europe** (Europe, Middle-East and Africa): +31 (0)320 277 070